

GLAND SELECTION CHART

CABLE SIZE	SINGLE CORE	2 CORE	3 CORE	4 CORE	5 CORE	7 CORE	10 CORE	12 CORE	19 CORE	27 CORE	37 CORE	48 CORE
1.5	-	16	16	20S.	20S.	20	20	25	25	32	32	32
2.5	-	20S.	20S.	20S.	20	20	25	25	25	32	40	40
4.0	-	20S.	20S	20	20	20	25	25	32	40	40	50
6.0	-	20	20	20	20							
10.0	-	20	20	25	25							
16.0	-	25	25	25	25							
25.0	-	25	32	32	32							
35.0	-	32	32	32	40							
50.0	20	25	32	32	50							
70.0	20	32	32	40	50							
95.0	25	32	40	50S								
120	25	40	50S	50								
150	32	40	50S	50								
185	32	50S	50	63S								
240	32	50	63S	63								
300	40	63S	63	75S								
400	40	63S	75S	75								

THE GLAND SIZES RECOMMENDED IN THIS CHART ARE FOR GUIDANCE ONLY AND ARE BASED ON CABLES MANUFACTURED TO BS 5467.

ACTUAL CABLE DIMENSIONS MUST BE USED FOR GLAND SELECTION.

